

Pride Colt Twin



Brand: Pride
Product Code: PC TWIN
Availability: 20
Weight: 61.00kg

Price: £2,939.00

Short Description

4mph pavement mobility scooter motability available.

Description

The sleek, sporty Colt Twin scooter delivers high performance operation and features unique dual front wheels to deliver the maximum stability of a four wheel scooter with the superb manoeuvrability of a three wheel scooter.

With a full complement of advanced standard features like exclusive low-profile tyres, wraparound delta tiller, and feather-touch the Colt Twin is a greater value than ever before.

Specification

- Weight Capacity 25 stone (158.75kg)
- Speed 4mph
- Ground Clearance: 2.25" (5.7 cm)
- Overall Length: 46" (116 cm)
- Overall Width: 22" (56 cm)
- Turning Radius: 42" (106.7 cm)
- Tyres: Front: 9" (22 cm)
- Rear: 10" (25 cm)

- Anti Tips – Rear
- Transport Weight: Total Weight (without batteries) 138lb (62.5Kg)
- Drive System: Rear wheel drive, sealed transaxle, 24 - volt DC motor
- Braking System: Electronic, regenerative, and electro-mechanical
- Battery requirements: 2 x 12 volt - 35 amp

[Read more about VAT Exemption.](#)

Product Gallery

