

Shoprider Perrero



Brand: Shoprider
Product Code: S888SLR
Availability: 20
Weight: 94.00kg

Price: £1,769.00

Short Description

6mph road/path friendly mobility scooter.

Description

Shoprider Perrero 6 mph mobility scooter:

- 1. Full suspension
- 2. Captain seat - Lockable four position swivel
- 3. Forward, back and height adjustable seat
- 4. High backrest - adjustable angle and fold down
- 5. Full light set / indicators
- 6. Horn and basket
- 7. High / low speed switch
- 8. Dual position armrest option designed to move with the seat or remain fixed on seat frame
- 9. Narrow profile wheels for easier steering
- 10. Adjustable tiller (five angle positions)
- 11. Height and width adjustable arms
- 12. Flip-up padded armrests
- 13. Anti-tippers
- 14. Wing mirror

Specification

- All round suspension
- Extra width seat
- 25 miles range
- 21 stone user weight
- Lights and indicators
- Battery Size - 2 x 12v 36ah
- Height - 121 cm (47.5)
- Length - 128 cm (50")
- Maximum Range - Up to 25 miles
- Maximum Safe Slope - 8 degrees
- Maximum Speed - 6 mph
- Maximum User Weight 136kgs (21st)
- Weight - 94kgs (207lbs)
- Width - 60 cm (24")

[Read more about VAT Exemption](#)

Product Gallery

