

Shanklin Mobile Footrest



Brand: Drive Medical
Product Code: 201
Availability: 15
Weight: 0.00kg

Price: £48.98

Short Description

Robust footrest enables users to raise and support their legs while sitting.

Description

Shanklin Mobile Footrest enables users to raise and support their legs while sitting

- Can be easily wiped clean
- Robust steel frame has a tough coating for durability
- Top size: 28cm x 50cm (11" x 19.7")
- Maximum weight capacity 100kg (16st)

Product Features

- Designed to enable users to raise and support their legs while sitting
- Comes with four 5cm (2") castors for ease of manoeuvring into position or out of the way for standing up from chair
- Horseshoe clips allow unit to be easily adjusted in both height and angle to accommodate chair height and pressure distribution requirements of legs/ankles
- Padded and contoured calf support allows comfortable and even distribution of pressure in legs

- Vinyl upholstery is fire retardant to IGM source 5 and can be easily wiped clean
- Robust steel frame has a tough coating for durability
- Top size: 28cm x 50cm (11" x 19.7")
- Height adjustable between 39cm - 54cm (15" - 21")
- Maximum weight capacity 100kg (16st)
- Please allow for 5-7 working days on all Welland products