

## Carisbrooke Static Footrest



**Brand:** Drive Medical  
**Product Code:** 200  
**Availability:** 15  
**Weight:** 0.00kg

**Price: £56.39**

### Short Description

Enables users to raise and support their legs while sitting.

### Description

Carisbrooke Static Footrest supports legs while sitting.

- Raises and supports legs while sitting
- Easily adjusted in both height and angle
- Padded and contoured calf support
- Robust steel frame has a tough coating

### Product Features

- Designed to enable users to raise and support their legs while sitting
- Horseshoe clips allow units to be easily adjusted in both height and angle to accommodate chair height and pressure distribution requirements of legs/ankles
- Padded and contoured calf support allows comfortable and even distribution of pressure in legs
- Vinyl upholstery is fire retardant to IGM source 5 and can be easily wiped clean
- Robust steel frame has a tough coating for longevity of use

- Top size: 28cm x 50cm (11" x 19.7")
- Comes with four 5cm (2") non-slip ferrules for secure positioning
- Maximum weight capacity 100kg (16st)
- Please allow for 5-7 working days on all Welland products