

## Breeze midi 4



**Brand:** TGA  
**Product Code:** BREEZEMIDI4  
**Availability:** 5  
**Weight:** 125.00kg

**Price: £3,848.00**

### **Short Description**

The new TGA Breeze Midi 4 brings to market a higher level of design features that now provide further comfort, reliability and safety.

### **Description**

TGA Breeze Midi 4 Building on the proven heritage of the popular larger Breeze models, this new Tga Breeze midi 4 version incorporates enhanced sharper looks, a disc braking system, fully orthopaedic seat and energy absorbing front and rear bumpers.

The tiller encompasses Delta 'Wig Wag' control levers and a gas sprung actuator which is easily adjustable for optimum ergonomic positioning for the owner. For fast status recognition, the easy-to-read LCD control panel displays all the relevant information the user requires in a simple format and safety is enhanced by a water-protected ignition switch.

If technical updates need to be uploaded or parameters changed, then the onboard computer can be accessed via the charging socket.

### **Specifications**

- Length 1380mm 54.33 inches

- Width 630mm 24.8 inches
- Turning radius 1470mm 58 inches
- Seat (Fully adjustable swivel seat)
- Total weight 125kg 276lbs
- Speed Variable up to 12km/h 8mph
- Approx range Standard battery 32km 20miles
- Battery size Standard battery 2 x 12v 50a/h
- Max carry capacity 130kg 20.4stone
- Max gradient 17%
- Ground clearance 90mm 3.5inches

[Read more about VAT Exemption](#)

## Product Gallery

