

Shoprider Granada



Brand: Shoprider
Product Code: S743
Availability: 5
Weight: 110.00kg

Price: £2,495.00

Short Description

Elegant and stylish road legal mobility scooter, max 25 miles range and user weight up to 25 stone with full suspension.

Description

Shoprider Granada is an elegant and stylish road legal mobility scooter.

Full Suspension

Features

- Max speed 8mph
- Range: 25 miles
- Max user weight: 25 stone
- Attractive removable modesty panel - more protection against elements of outdoors
- LED dashboard for greater clarity
- Full LED lights set - front and rear
- Dual Mode wig-wag lever - greater options on how to move the scooter
- Automatic half speed sensor on full lock turns in 8mph setting for additional safety
- Rear disc brake

- Non-marking pneumatic tyres
- Comfortable adjustable height, reclining captain seat
- Flip up, width and height adjustable armrests
- Infinite angle tiller adjustment
- Front shopping basket with carry handle

Specifications

- Overall length 130cm
- Overall width 66cm
- Overall height 124cm
- Seat width 46cm
- Seat depth 41cm
- Seat height 60-72cm
- Backrest height from seat 55cm
- Seat to footrest/floorpan height 40-52cm
- Safe working slope 8°
- Maximum user weight 160kg/25stone
- Battery size 12v x 2 50ah
- Range 25 miles
- Weight(with battery pack, seat, basket) 110kg
- Maximum speed 8mph
- Wheel/tyre size 30cm
- Turning radius 126cm
- Battery charger 5 amp off board
- Ground clearance 13cm

 [Read more about VAT Exemption...](#)