

## Traveso



**Brand:** Roma Medical  
**Product Code:** 889 XLSBN CAN/SILVER  
**Availability:** 7  
**Weight:** 195.00kg

**Price: £8,495.00**

### **Short Description**

This Shoprider Traveso, the only scooter on the market like it, looks like a scooter which could be classed as a car, and that's what it is.

### **Description**

The Shoprider Traveso looks like a scooter which could be classed as a car, and that's what it is, comes with sleek body work and stylish curves this is the only scooter on the market like it.

Durable rustproof cover to keep you dry. The windows can be opened, and the doors can even be removed to give you ventilation as well as safety and protection.

The windscreen wipers and rear view mirrors are also fitted as standard, and there is a handbrake for your additional safety. This scooter is also lockable so it can be left safely.

Incorporates a weatherproof environment without being exposed to the elements. The features, comfort and performance are unrivalled.

A Hi / Lo speed switch is fitted to reduce the speed to 4mph so that the scooter can be used legally on the pavement as well as on the road.

The Delta handlebar means that the scooter can be steered and controlled using either the left or right hand. The adjustable Captain seat provides a superb level of comfort, and allows you to find your perfect driving position easily.

All these features are combined with advanced steering controls and suspension.

Now you can go out or return home in the worst of weathers, in comfort and safety.

## Features

- Maximum range of 20 miles
- Maximum user weight of 35 stone / 222Kg
- Full light set
- Rustproof all weather lockable cap
- Delta tiller controls
- Swivelling Captain seat with headrest
- Full suspension
- 13" / 33cm wide profile wheels
- Handbrake
- Doors on both sides
- Wiper and rear view mirror as standard
- Genuine automotive windshield for climate safety

## Specifications

- Adjustable speed
- Height 162 cm
- 160 cm
- Maximum range 20 miles
- Maximum user weight 222 kg
- Number of wheels 4
- Road legal
- Safe climbing angle 12°
- Weight with batteries 195 kg
- Width 73 cm

 [Read more about VAT Exemption...](#)