

Breeze midi 3



Brand: TGA
Product Code: BREEZEMIDI3
Availability: 5
Weight: 110.00kg

Price: £3,848.00

Short Description

The new compact version of the Breeze 3 features updated looks, disc brakes, orthopaedic seat and new energy absorbing bumpers.

Description

The updated disc brakes of Breeze midi 3 scooter are more responsive and incorporate an emergency stopping system if the user loses control of the scooter.

A set of 50ah batteries will allow the user to enjoy journeys over 34miles to a full charge, with a ground clearance of 3.5in and sturdy construction it allows you to tackle terrain that is not the smoothest.

Similar to its 4-wheel sister, the easily manoeuvrable Breeze Midi 3 encompasses full front and rear suspension which delivers a smooth, sturdy ride across both even and rough surfaces whether travelling at slow speeds or its top limit of approximately 12.8km/h (8mph).

Powered by 2 x 50A batteries, the Breeze Midi 3 has a range of 50-55km (31.1-34.2miles), can drive safely up a 9.5% incline and an impressive ground clearance of 90mm (3.54") which is ideal for driving off road in a more rural location.

Specifications

- Length 1330mm 51inches
- Width 675mm 25inches
- Turning radius 990mm 39inches
- Seat (Fully adjustable swivel seat)
- Total weight 110kg 243lbs
- Speed Variable up to: 12km/h 8mph
- Approx range Standard battery 32km 20miles
- Battery size Standard battery 2 x 12v 50a/h
- Max carry capacity 130kg 20.4stone
- Max gradient 17%
- Ground clearance 90mm 3.5inches

[Read more about VAT Exemption](#)

Product Gallery

