

## Deluxe Handy Reacher



**Brand:** Aidapt  
**Product Code:** VM901F/H/K/M  
**Availability:** 30  
**Weight:** 0.00kg

**Price: £8.18**

### Short Description

Deluxe Handy Reacher features an improved gripping claw and comfortable handle/trigger.

### Description

The Deluxe Handy Reacher is part of the Aidapt range of personal aids and features an improved gripping claw and comfortable handle/trigger. Handy Reachers can help to reduce excessive pressure on your back, as sometimes experienced when reaching or bending over. A magnet is attached to the tip for picking up smaller items.

- Featuring a 360 degrees rotatable head
- Magnetized tip for retrieving small metallic items
- Lightweight and durable construction for ease of use

	<b>VM901F</b>	<b>VM901H</b>	<b>VM901K</b>	<b>VM901M</b>
<b>Depth (mm)</b>	40	40	40	40
<b>Height (mm)</b>	610	760	610	890
<b>Width (mm)</b>	125	125	125	125
<b>Product Dimensions (mm)</b>	610x125x40	760x125x40	810x125x40	890x125x40

<b>Net weight (kg)</b>	0.18	0.18	0.18	0.18
<b>Colour</b>	Black and chrome	Black and chrome	Black and chrome	Black and chrome
<b>Magnet</b>	yes	yes	yes	
<b>Size</b>	Length: 600 mm (24 inch)	Length: 750 mm (30 inch)	Length: 800 mm (32 inch)	Length: 875 mm (35 inch)

[Read more about VAT Exemption...](#)