

Colt Sport



Brand: Pride
Product Code: PCSPORT
Availability: 6
Weight: 68.00kg

Price: £3,549.00

Short Description

The all-new compact Colt Sport Scooter features a sleek, sporty design and delivers high-performance operation.

Description

Pride Colt Sport standard features include a handbrake, front and rear suspension, feather-touch disassembly, a wraparound delta tiller, high-visibility automatic brake lights and an LED kerb light.

With a full complement of advanced standard features including, rugged suspension, high-visibility automatic front and rear lighting system, feather-touch disassembly, LED kerb light, and wraparound, delta tiller, the Colt Deluxe electric disabled scooter is a greater value than ever before.

- Max Range: 26.8 miles
- Max Weight Capacity: 25 stone
- Ground Clearance: 8.25cm (3.25")
- Turning Radius: 133.35cm (52.5")
- Overall Length: 119.38cm (47")
- Overall Width: 66.37cm (26.13")
- Front Tyres: 9 cm x 25.4 cm (3.5" x 10")
- Rear Tyres: 10cm x 25.4cm (4" x 10")

- Suspension: Front and rear
- Maximum obstacle climbing ability: 5cm (2")
- Total Weight without batteries: 68kg (150lbs)
- Heaviest Pieces when disassembled (front section): 29.48kg (65lbs)
- Seat Weight: 19.05kg (42lbs)
- Seat Dimensions: 44.45cm wide x 43.18cm deep
- Drive system: rear-wheel drive, sealed transaxle, 24 volt DC motor
- Dual braking system: regenerative and electromechanical
- Battery Requirements: 2 x 12 Volt, deep cycle
- Battery Size: 45ah
- Battery Weight: 14.74kg or 32.5lbs
- Battery Charger: 5amp off board
- Available in choice of Red or Blue colours

[Read more about VAT Exemption](#)

Product Gallery



