

Shoprider Cadiz



Brand: Roma Medical
Product Code: S889SL
Availability: 7
Weight: 100.00kg

Price: £2,495.00

Short Description

Powerful 8mph road legal mobility scooter with full Suspension and full lighting package.

Description

Shoprider Cadiz is a powerful 8mph road legal mobility scooter with full suspension, full lighting package, max 20 miles and user weight up to 25 stone.

Features

- Captain seat- lockable four position swivel
- Forward, back and height adjustable seat
- Width adjustable arms
- Full light set
- High/low speed switch
- Front and rear suspension
- Adjustable angle tiller with delta handle
- Large pneumatic tyres for extra comfort
- Available in red and black (additional cost for black)

Specifications

- Overall length 1300mm/51"
- Overall width 640mm/25"
- Overall height 1120mm/44"
- Seat width 460mm/18"
- Seat height (to ground) 630-705mm/25"-28"
- Armrest height 210mm/8"
- Backrest height from seat (including headrest) 500mm-650mm/20"-26"
- Seat to footrest/floor height 430mm-505mm/17"-20"
- Safe working slope 8°/14%
- Maximum user weight 160kg/25st
- Battery size 12v x 2 36Ah (45Ah optional £75)
- Range 20 miles (32km)
- Weight(with battery pack, seat, basket) 100kg/220 lbs
- Maximum speed 8mph
- Wheel/tyre size 330mm/13"
- Turning radius 1500mm/59"
- Battery charger 3 amp off board
- Ground clearance*** 150mm/6"

[Read more about VAT Exemption.](#)

Product Gallery



