

Perching stool



Brand: Aidapt
Product Code:
VG862/VG863/VG864/VG865
Availability: 20
Weight: 0.00kg

Price: £56.10

Short Description

Enables people who are unable to stand for long periods of time to work in the kitchen, or wash in the shower / at the bathroom basin.

Description

The Astral Perching Stool is the first mass-produced perching stool range to include a unique, seamless polyurethane seat. Developed for comfort and longevity, the contoured lines of both the seat and back are not only attractive to the eye, but also offer excellent ergonomic support. The height adjustable frame has side mounted stress bars and a large footprint for safety and stability. The Astral Perching Stool is available with a range of optional fittings.

- Large footprint for stability
- Ideal for those who struggle to stand for long periods
- Maximum user weight limit 159kg (25 stone)
- Height adjusting legs to suit most users
- Padded wipe clean PU seat for extra comfort

	VG862	VG863	VG864	VG865
Depth (mm)	580	580	585	585
Height (mm)	500-650	-	500-680	790-940
Width (mm)	510	510	410	410

Product Dimensions (mm)	500x510x580		410x585x500	790x410x585
Net weight (kg)	2.8	2.8	2.8	2.8
Back Rest (mm)	-	-	-	350x150x25
Colour	White and Blue	White and Blue	White and Blue	White and Blue
Configuration	Plain Perching Stool	With Arms Only	With Arms and Plain Back	With Arms and Padded Back
Footprint (WxD) Min/Max (mm)	510x335	510x335	410-460x585-650	410-460x585-650
Footprint Depth (mm)	335	335	585-650	585-650
Height of Arms-from Seat (mm)		-	170	170
Maximum User Weight (kg)	159	159	159	159
Maximum User Weight (stones)	25	25	25	25
Overall Height (mm)		790-940	500-680	790-940
Padded Seat (WxD) (mm)	360x290	360x290	350x300	350x300
Seat Height (mm)	500-650	500-650	30	500-650
Seat Width (mm)	360	360	350	350
Width between Arms (mm)	490	490	480	480

 [Read more about VAT Exemption...](#)

Product Gallery

