

## Furniture Raisers Grip on



**Brand:** Aidapt  
**Product Code:** VG823D/VG823F/VG823H  
**Availability:** 25  
**Weight:** 0.00kg

**Price: £47.97**

### Short Description

Suitable for furniture on castors or legs, raisers 'Grip on' using the unique 'donut'.

### Description

Grip-on Chair Raisers, suitable for use with furniture on castors or legs, these raisers grip on to the furniture using the unique 'donut'. This greatly increases safety as the furniture raiser is attached to the furniture rather than the furniture simply sitting on the raiser, as with some other models. Raisers are supplied in a set of 4. To raise 25 mm (1") one set is required, 2 sets are required to raise 50 mm (2") and 3 sets are required to raise 75 mm (3"). Available to fit several widths of castors/legs. Maximum user weight 500 kg (78 stone) (including furniture).

- Raises furniture to aid bending when getting in and out of beds and chairs
- Suitable for use with furniture on castors or legs
- Supplied as a set of 4
- Maximum user weight limit 500kg (78 stones) including furniture

	<b>VG823D</b>	<b>VG823F</b>	<b>VG823H</b>
<b>Depth (mm)</b>	92	100	92
<b>Height (mm)</b>	120	120	120

<b>Width (mm)</b>	92	100	92
<b>Product</b>	120x92x92	120x100x100	
<b>Dimensions (mm)</b>			
<b>Net weight (kg)</b>	0.275	0.275	0.275
<b>Maximum User</b>	500	500	500
<b>Weight (kg)</b>			
<b>Maximum User</b>	78	78	78
<b>Weight (stones)</b>			
<b>Size</b>	Max. Diameter: 32 mm (1 1/4")	Max. Diameter: 45 mm (1 3/4")	Max. Diameter: 58 mm (2 1/4")

[Read more about VAT Exemption](#)

## Product Gallery

