

## The Lars Luxurious Leather



**Brand:** Drive Medical  
**Product Code:** ALC  
**Availability:** 120  
**Weight:** 0.00kg

**Price: £1,219.00**

### **Short Description**

The Lars Riser Recliner is upholstered in modern and luxurious leather with dual motor mechanism.

### **Description**

The dual motor will allow you to operate and move the leg rest and the back rest independently.

- Stylish and modern leather design
- Excellent snooze position
- It will rise and tilt users to their feet
- Quality assured

### **Product Features:**

- Stylish and modern leather design to add to any decor
- Two motors enable leg rest to be elevated / lowered and backrest to be reclined independently to give more flexible support options
- Back rest can be positioned upright while legrest is elevated, for example to watch TV
- It will rise and tilt users to their feet
- Excellent snooze position for supreme comfort and relaxation

- Quiet motors provide a smooth and discreet operation
- Ergonomically contoured seat and back and padded upholstery for added comfort and support
- Convenient side pocket for easy storage
- Safety key is built into hand controller to prevent accidental movement and entrapment
- Quality assured, meeting all relevant fire retardancy standards
- In built fail safe battery back-up system which enables the chair to be raised and lowered in the event of a power failure\*
- Dual motor comes with a 5 button handset
- Features castors to aid movement across all surfaces
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems
- Available in luxurious lather with 4 colour options – Black, Burgundy, Brown or Cream

**Specifications: Dimensions:** (W x D x H): 98cm (38.5") x 95cm (37.5") x 111cm (43.5")

**Seat:** (W x D x H): 50cm (20") x 48cm (19") x 50cm (20")

**Weight Capacity:** 135kg (21st)

[Read more about VAT Exemption](#)

## Product Gallery



