

Dual Motor Fabric Chair



Brand: Drive Medical
Product Code: CLR32
Availability: 60
Weight: 0.00kg

Price: £745.00

Short Description

The Dual Motor features a mechanism that allows independent operation of the leg rest and back rest.

Description

The Dual Motor Rise & Recliner features a mechanism that allows independent operation of the leg rest and back rest. This provides countless recline options for finding the most comfortable position.

- Dual Motor for flexible support options
- Quiet and simple operation
- Ergonomically contoured
- Convenient side pocket

Product Features:

- Two motors enable legrest to be elevated / lowered and backrest to be reclined independently to give more flexible support options
- Full length chaise
- Back rest can be positioned upright while leg rest is elevated, e.g. example to watch TV
- Three tier waterfall back design provides extra comfort and support

- It will rise and tilt users to their feet
- Quiet motors provide a smooth and discreet operation
- Ergonomically contoured seat and back and padded upholstery for added comfort and support
- One convenient side pocket for easy storage
- Safety key is built into hand controller to prevent accidental movement and entrapment
- Quality assured, meeting all relevant fire retardant standards
- Comes with a simple to operate 4 button handset
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems

Specifications:

Dimensions: (W x D x H): 72cm (28") x 85cm (33") x 103cm (40.5")

Seat: (W x D x H): 53cm (21") x 54cm (21") x 51cm (20")

Weight Capacity: 140kg (22st)

[Read more about VAT Exemption](#)

Product Gallery



