

## Three tier PU Riser (Dual)



**Brand:** Drive Medical  
**Product Code:** CLR33  
**Availability:** 60  
**Weight:** 0.00kg

**Price: £735.00**

### **Short Description**

The Dual Motor Rise & Recliner features a mechanism that allows independent operation of the leg rest and back rest. This provides countless recline options for finding the most comfortable position.

### **Description**

The two motors of the Three Tier PU Riser enable leg rest to be elevated / lowered and back rest to be reclined independently to give more flexible support options

It is fully padded for increased comfort, including a three-tier waterfall back.

- Two motors for flexible support options
- Easy to clean PU material
- Quiet and simple operation
- Ergonomically contoured chair

### **Product Features:**

- Two motors enable leg rest to be elevated / lowered and backrest to be reclined independently to give more flexible support options
- Easy to clean soft feel PU material full length chaise
- Backrest can be positioned upright while leg rest is elevated, e.g.

example to watch TV

- Three tier waterfall back design provides extra comfort and support
- It will rise and tilt users to their feet
- Quiet motors provide a smooth and discreet operation
- Ergonomically contoured seat and back and padded upholstery for added comfort and support
- One convenient side pockets for easy storage
- Safety key is built into hand controller to prevent accidental movement and entrapment
- Quality assured, meeting all relevant fire retardant standards
- Comes with a simple to operate 4 button handset
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems

### Specifications:

**Dimensions:** (W x D x H): 72cm (28") x 85cm (33") x 103cm (40.5")

**Seat:** (W x D x H): 53cm (21") x 54cm (21") x 51cm (20")

**Weight Capacity:** 140kg (22st)

**COLOURS:** Brown, Cream, Black

[Read more about VAT Exemption](#)

### Product Gallery



