

## Scooter Canopy Fold Away



**Brand:** 1st Choice Mobility  
**Product Code:** SCOOTERPAC  
**Availability:** 6  
**Weight:** 0.00kg

**Price: £545.00**

### **Short Description**

The world's first fold-away mobility scooter canopy.

### **Description**

No long waits for a canopy to be made, and with the added benefit of it fitting virtually every scooter-instantly. From the largest 8mph scooters right through to the smallest car-transportable the ScooterPac clever design means one size fits all.

The ground breaking patent-pending technology provides protections from the sun and lets you soak up the sun (when its out!) If the scooter is exchanged then the ScooterPac is easily swapped to the new scooter - with no waiting around! Aluminium construction for long term durability and no chance of rust! And state of the art cover material too!

Fitting is simple too! Simply slide the bracket into the accessory bar on the bottom of the scooter seat, and attach the adjustable arms to the scooter armrests.

The ScooterPac canopy is unique, comprising a light yet strong and long-lasting aluminium frame, hard wearing waterproof canvas roof, clear polypropylene side windows and a hard, clear Perspex front window.

Suitable for all types of mobility scooter.

[Read more about VAT Exemption](#)

## Product Gallery

