

## Walking Stick Ice Grip



**Brand:** Drive Medical  
**Product Code:** WSA006RT  
**Availability:** 30  
**Weight:** 0.00kg

**Price: £9.99**

### Short Description

Walking Stick Ice Grip provides grip in snow or icy weather conditions.

### Description

#### Specifications:

- Five pronged Ice Grip assists users when using their walking aids in snow or icy weather conditions
- Suitable for most walking sticks and crutches with tubing / ferrule sizes between 19mm and 22mm
- Ice grip can be 'flipped up' when not in use into a locked position to allow walking stick to be used as normal
- Supplied with 2 pairs of screws; shorter screws to install the ice grip on narrow tubing and longer screws for wider tubing

 [Read more about VAT Exemption](#)

### Product Gallery

