

## Charger In car



**Brand:** STRIDENT Batteries  
**Product Code:** CHHP24-1.5CAR  
**Availability:** 10  
**Weight:** 0.00kg

**Price: £75.00**

### Short Description

In car charger charges 7ah – 15ah batteries.

### Description

For convenient use from the car charging point, this 24V 1.5Amp High Power in-Car Mobility Scooter Charger has been designed for use with Mobility Scooters and Powerchairs with 7Ah to 15Ah batteries.

It employs the Strident three stage charging technique that greatly reduces premature battery fatigue and dramatically improves performance for complete reassurance and peace of mind for the user.

The charger is quick and easy to set up. It conditions and maintains the scooter's battery for premium functionality and reliability.

### Product Specification

- One 24V 1.5Amp High Power In-Car Mobility Scooter Charger
- Suitable for charging 7Ah - 15Ah batteries
- Compatible with most mobility scooters
- Intelligent three stage charging technique
- Ensures optimum battery performance

- Internal cooling fan
- Strong, durable metal casing
- Provides years of reliable service
- Quick and easy to use
- Case length: 132mm (5")
- Case width: 92mm (3.5")
- Case height: 53mm (2")

 [Read more about VAT Exemption...](#)