

Ashby Raised Toilet Seat



Brand: Aidapt
Product Code: VR207
Availability: 25
Weight: 0.00kg

Price: £42.39

Short Description

If you struggle when lowering yourself onto the toilet the raised toilet seat (2", 4" and 6") reduces stress on your knees and hips...

Description

The Ashby Super Contract Raised Toilet Seat is a hygienic one-piece moulded seat. The Ashby Raised Toilet Seat is an ideal product for trade/volume users; the design lends itself to easy cleaning and reissue. Additionally the Ashby Toilet Seat has an integral anti-bacterial finish to reduce the potential of cross infection. The Ashby Toilet Seat has a front mounted and two rear adjustable brackets; each bracket is made from steel and coated with a rubber sleeve to maximise grip and stability.

Specifications

Code	VR207	VR208	VR209
Aperture Depth (mm)	250	250	250
Aperture Width (mm)	220	220	220
Colour	White	White	White
Depth (mm)	250	250	250

Footprint (WxD)	220x250	220x250	220x250
Min/Max (mm)			
Footprint Depth (mm)	250	250	250
Footprint Width (mm)	220	220	220
Height (mm)	51	102	152
Maximum User Weight (kg)	190	190	190
Maximum User Weight (stones)	30	30	30
Net weight (kg)	0.95	1	2
Overall Width (mm)	220	220	220
Packing	Poly Bag		
Product Dimensions (mm)	51x220x250	102x220x250	152x220x250
Seat Depth (mm)	250	250	250
Seat Height (mm)	51	102	152
Seat Width (mm)	220	220	220
Weight Limit (kg)	190	190	190
Width (mm)	220	220	220
Width of Front Cutaway (mm)	78	78	78
Size	Seat Height: 50 mm (2")	Seat Height: 100 mm (4")	Seat Height: 150 mm (6")

[Read more about VAT Exemption](#)

