

Suction grab bar rail



Brand: Aidapt
Product Code: VY460
Availability: 30
Weight: 0.00kg

Price: £13.27

Short Description

The Quick Suction Grab Rail is a compact support handle and offers greater independence by providing stability when getting in and out of the bath or shower, or for general use around the home.

Description

It is ideal for those who need a little extra mobility support. You can install the product in seconds, with no need for messy drilling or DIY, by simply placing the grip onto any smooth, non-porous surface (such as tiles) and pushing the safety locking levers down. For total confidence, the device will indicate it is securely locked in place with a green indicator on the suction handle. It is a great product to have when travelling and staying with friends or relatives or in hotels; providing a non-permanent support solution, without damaging the surface below. NB. This product is designed to assist your natural motion, as when getting in and out of a shower. It is not designed to hold or support your full body weight.

Specifications

- Designed for use on tiled walls or other non-porous surfaces
- Ideal for use when travelling/staying with relatives
- Can be removed in seconds for easy cleaning
- No messy drilling or screwing required to fit

- Can help prevent dangerous slips and falls
- Powerful, secure locking suction pads with safety indicator
- Colour: White/Blue
- Depth (mm): 75
- Depth to Wall (mm): 80
- Height (mm): 292
- Length (Extended) (mm): 295
- Net weight (kg): 0.35
- Product Dimensions (mm): 292x110x75
- Width (mm): 110

SAFETY NOTICE

Suction Grip Grab Rails are designed for user stability only. They are not intended to support full body weight. Ratings shown are for maximum load.

Product Gallery



