

President Grab Bar



Brand: Aidapt
Product Code: VY425/VY426/VY427
Availability: 30
Weight: 0.00kg

Price: £7.43

Short Description

The President Grab Bar from Aidapt is an attractive, durable grab bar constructed from UPVC and finished in white; the 32 mm tube is ribbed for improved grip.

Description

The President Grab Bar is an attractive, durable grab bar ideal for volume use, with a 32 mm diameter bar. It is constructed from UPVC and the tube is ribbed for improved grip. The flanges are multi-fixing giving the choice of five different fixing points to suit all installation needs. Each fixing point is reinforced and the flange has a cover cap as standard. This model is finished in white and is available in several lengths.

- Ribbed grip for improved safety
- Reinforced fixing points
- No earth bonding required
- Cover caps as standard
- Ideal for bathroom walls for extra security when getting in and out of the bath
- Product Dimensions (mm): 292x110x75
- Width (mm): 110
- **Wall fixings not supplied**

	VY425	VY426	VY427
Depth (mm)	100	100	-
Height (mm)	-	-	670
Width (mm)	80	80	90
Product Dimensions (mm)	380x80x100	530x80x100	670x90x75
Net weight (kg)	0.31	0.31	0.31
Colour	White	White	White
Depth to Wall (mm)	60	60	60
Size	Length: 300 mm (12")	Length: 450 mm (18")	Length: 600 mm (24")
Tube Diameter (mm)	32	32	32

SAFETY NOTICE

Suction Grip Grab Rails are designed for user stability only. They are not intended to support full body weight. Ratings shown are for maximum load.

[Read more about VAT Exemption...](#)