

Deluxe Pressure Relief Orthopaedic Coccyx Cushion

Brand: Aidapt
Product Code: VA126WA
Availability: 20
Weight: 0.00kg



Price: £49.49

Short Description

Deluxe Pressure Relief Orthopaedic cushion is constructed from a 3 cm thick high density foam topped with a generous 4.5 cm eggcrate (convoluted) memory foam topper for improved comfort and air circulation

Description

Bring a little comfort into your life with the Deluxe Pressure Relief Orthopaedic Cushion. These cushions are designed to provide extra support and comfort in cars, homes or wheelchairs. Featuring a black cotton twill cover with two pairs of straps to secure the cushion to a wheelchair or dining chairs, for example. The Deluxe Pressure Relief Orthopaedic Cushion is constructed from a 3 cm thick high density foam topped with a generous 4.5 cm memory foam topper. The convoluted (eggcrate) memory foam helps to stimulate circulation while reducing pressure and discomfort. It also allows more air circulation, keeping you cool when using the cushion for longer periods.

- Adds comfort for regular wheelchair users
- Removable washable cover
- Can be secured to wheelchairs using integrated tie handles
- Ideal for providing extra support in cars and around the home

- Eggcrate (convoluted) memory foam design for improved air circulation
- Height (mm): 85
- Width (mm): 430
- Product Dimensions (mm): 85x430x430
- Net weight (kg): 0.75
- Maximum User Weight (kg): 100
- Maximum User Weight (stones): 15

[Read more about VAT Exemption](#)

Product Gallery

