

Vinyl Wheelchair Cushion with Memory Foam

Brand: Aidapt
Product Code: VA122SA
Availability: 20
Weight: 0.40kg



Price: £16.34

Short Description

Vinyl Wheelchair Cushion with Memory Foam is designed to make you more comfortable.

Description

Bring a little comfort into your life with the Vinyl Cushion for Wheelchairs. These wheelchair cushions are designed to make your ride more comfortable; the inner foam has a generous memory foam topper for a touch of luxury. The wipe-clean, black vinyl cover with zip closure features ties at the two back corners to quickly and easily secure the cushion in place on wheelchairs or other chairs. Available in various sizes to fit most types of wheelchair, the cushion can also provide extra support and comfort in cars and homes.

- Ideal for use in the car or around the home
- Integrated ties secure the cushion to the wheelchair frame
- Generous memory foam top layer adds a little luxury
- Super-soft vinyl cover with zip closure is easy to clean
- Enjoy extra comfort on most types and sizes of wheelchair

	VA122SA	VA123SA	VA124SA	VA125SA
Depth (mm)	50	50	50	50
Height (mm)	406	457	406	460

Width (mm)	406	406	406	410
Product Dimensions (mm)	50x406x406	457x406x50	406x406x75	460x410x78
Net weight (kg)	0.4	0.4	0.4	0.4
Colour	Black	Black	Black	Black
Main body material	Memory/PU Foam/Vinyl	Memory/PU Foam/Vinyl	Memory/PU Foam/Vinyl	Memory/PU Foam/Vinyl
Size	406x406x50 mm (16x16x2 inches)	457x406x50 mm (18x16x2 inches)	406x406x75 mm (16x16x3 inches)	457x406x75 mm (18x16x3 inches)
	VA122SA	VA123SA	VA124SA	VA125SA
Product EAN	5021196699651	5021196699668	5021196699675	5021196699699
Packaged dimensions (HxWxD mm)	410x410x52	460x410x52	410x410x78	460x410x78
Display Footprint (mm)	410x410x52	460x410x52	410x410x78	460x410x78
Transport dimensions (mm)	410x410x52	460x410x52	410x410x78	460x410x78
Gross Weight (Packaged) (kg)	0.8	0.8	0.8	0.8
Packaging	Colour Box	Colour Box	Colour Box	Colour Box
Item priced per	1	1	1	1
Inner box quantity	0	0	0	0
Outer box (HxWxD mm)	600x430x440	610x490x450	810x430x440	850x490x450
Outer box quantity	9	10	10	10
Gross weight outer carton (kg)	8	8	8	8
Amazon Asin	B00CIW4IZG	B00CIW4KP4	B00CIW4LEY	B00CIW4MCA

[Read more about VAT Exemption](#)

Product Gallery

