

## Deluxe Pressure Relief Orthopaedic Coccyx Cushion



**Brand:** Aidapt  
**Product Code:** VM974AB  
**Availability:** 19  
**Weight:** 0.00kg

**Price: £28.04**

### Short Description

Deluxe Pressure Relief Orthopaedic Coccyx Cushion is constructed from a 3 cm thick high density foam topped with a generous 4.5 cm eggcrate (convoluted) memory foam topper for improved comfort and air circulation.

### Description

Bring a little comfort into your life with the Deluxe Pressure Relief Orthopaedic Coccyx Cushion. These wheelchair cushions are designed to make your ride more comfortable and can also provide extra support and comfort in cars and homes. Featuring a black cotton twill cover with two pairs of straps to secure the cushion to a wheelchair or dining chairs, for example. The Deluxe Pressure Relief Orthopaedic Coccyx Cushion is constructed from a 3 cm thick high density foam topped with a generous 4.5 cm memory foam topper for improved comfort and air circulation. The convoluted (eggcrate) memory foam helps to stimulate circulation while reducing pressure and discomfort. It also allows more air circulation, keeping you cool when using the cushion for longer periods.

- Eggcrate (convoluted) memory foam design for improved air circulation
- Designed to fit most types of wheelchair
- Provides extra support in cars, around the home or on a wheelchair

- Depth (mm): 70
- Height (mm): 430
- Width (mm): 430
- Product Dimensions (mm): 430x430x70
- Net weight (kg): 0.525
- Aperture Depth (mm): 100
- Aperture Width (mm): 100
- Colour: BLACK
- Maximum User Weight (kg): 100
- Maximum User Weight (stones): 15

[Read more about VAT Exemption](#)

## Product Gallery

