

Sock helper



Brand: Aidapt
Product Code: VM947A
Availability: 30
Weight: 0.00kg

Price: £14.85

Short Description

Helps people to put on socks or stockings without bending down to the feet.

Description

The Rigid Sock & Stocking Helper makes putting on socks and tights easy. Sock Helpers are used by simply stretching the garment over the inner metal frame and then the foot can easily be inserted. Use the integral handles of the sock and panty hose aid help to then slide the foot into the sock/stocking. This helper is especially designed to help to reduce the need to bend and is ideal for users recovering from surgery and those with limited mobility.

- A great solution for those who struggle to bend when putting socks/stockings on.
 - Can be used to put on the majority of socks, stockings and tights.
 - Lightweight design
 - Easy to use
 - Helps to reduce the need to bend when dressing.
-
- Depth (mm): 170
 - Height (mm): 350
 - Width (mm): 205
 - Product Dimensions (mm): 350x205x170

- Net weight (kg): 0.75
- Warranty Period: 12 Months
- Colour: White

[Read more about VAT Exemption...](#)

Product Gallery

