

Button hook and zipper



Brand: Aidapt
Product Code: VM941
Availability: 30
Weight: 0.00kg

Price: £6.20

Short Description

Button hook and zipper enables people with reduced hand function to pull up zips and do up buttons independently.

Description

A handy device for securing buttons and fastening zips, the Plastic Handle Button Hook/Zipper has a loop hook at one end for buttons and a C hook at the other for zips. It is ideal for people who have limited mobility or the use of only one hand. It is available in red, black and blue.

- Dressing aid that helps users with arthritis or a weakened grip to fasten buttons and zippers.
- Can be used with the majority of buttons and zips.
- Lightweight and attractive design is ideal for travel.
- Available in three colours.
- Easy to use.

	VM941A	VM941B	VM941C
Depth (mm)	10	10	10
Height (mm)	200	200	200
Width (mm)	30	30	30
Product	200x30x10	200x30x10	200x30x10

Dimensions (mm)

Net weight (kg)

0.42

0.42

0.42

Colour

BLACK

RED

BLUE

[Read more about VAT Exemption](#)

Product Gallery

