

Pride Apex Lite



Brand: Pride
Product Code: APEXLIT
Availability: 20
Weight: 36.90kg

Price: £1,469.00

Short Description

Modern, Lightweight and agile the Apex Lite from Pride Mobility.

Description

The Apex Lite features 8" Alloy style wheels, detachable frame design and a comfortable folding seat aiding for a smoother more enjoyable journey. The Apex Lite is the perfect travel partner, helping you get where you need to go.

Specification

Battery Weight (each)	9.07 kg
Weight Capacity	21 Stone 6 lbs (136 kg)
Maximum Speed	4 mph (6 km/h)
Ground Clearance	3.4 cm (1.3")
Overall Length	102.2 cm (40.2")
Overall Width	49 cm (19.2")
Turning Radius	138.5 cm (54.4 cm)
Standard Controls	n/a
Heaviest Piece (when disassembled)	15 kg Front Section
Total Weight (without batteries)	36.9 kg
Suspension	None
Battery Requirements	Two 12-volt, deep cycle (standard) Size:

Seating
Tyres

2 x 12 AH
Foldable, Compact.
Front: 20 cm Rear: 20 cm

[Read more about VAT Exemption](#)

Product Gallery

