

## Pride Go-Go Elite Traveller LX



**Brand:** Pride  
**Product Code:** PGG LX  
**Availability:** 20  
**Weight:** 35.60kg

**Price: £1,679.00**

### Short Description

Go-Go Travel Scooters LX are the first choice in travel mobility for many thousands of customers the world over. Truly transportable and passionate loved, what can the Go Go Range offer you?

### Description

The Go-Go Elite Traveller® LX delivers luxury and convenient features like a wraparound delta tiller with high level charging port, LED kerb light and large storage basket together with a compact, easy to transport design.

Sleek styling complements the Go-Go Elite Traveller LX's abundant features making it an exceptional travel mobility choice.

### Model Specification

- Weight Capacity: 19 Stone 9 lbs
- Transport Weight: Total Weight (without battery pack) 35.6 kg (78.5 lbs)
- Heaviest Part ( front section) 15.2 kg (33.5 lbs)
- Ground Clearance: 2.50" (6.35 cm)
- Drive System: 24-volt DC motor, sealed transaxle, rear-wheel drive
- Turning Radius: 54.55" (138.5 cm)
- Dual Braking System: Regenerative and electromechanical

- Overall Length: 39.5" (100.33 cm)
- Standard Seating: Foldable, compact, vinyl-covered molded plastic
- Overall Width: 19.5" (49.5 cm)
- Battery Charger: Off-board, 2-amp
- Tyres: Front: 2" x 7" (5 cm x 17.78 cm)
- Back: (2.50" x 8") (6.35 cm x 20.32 cm)
- Anti Tips - Rear
- Battery Requirements: 2 x 12 volt - 12 amp batteries in one battery pack.

[Read more about VAT Exemption](#)

## Product Gallery



