

Pride Go-Go Elite Traveller Sport



Brand: Pride
Product Code: GOGOSPORT
Availability: 20
Weight: 44.00kg

Price: £1,949.00

Short Description

Go-Go Elite Traveller Sport is the ultimate high-performance scooter with added larger deck space.

Description

The Go-Go Elite Traveller Sport with added larger deck space, features an innovative design with a full complement of unique standard features. A charger port conveniently located in the tiller, front and rear LED lighting, 9-inch wheels at front and rear and an increased weight capacity of 23 Stone 3 lbs. are just a few features that make the Go-Go Elite Traveller Sport an exceptional scooter.

Specifications:

Specification Battery Weight (each)	13 kg (29 lbs.)
Weight Capacity	23 Stone 3 lbs (147 kg)
Maximum Speed	4 mph (6 km/h)
Ground Clearance	3.8 cm (1.5") at motor
Overall Length	109.2 cm (43")
Overall Width	54 cm (21.25")
Turning Radius	128.3 cm (50.5")
Standard Controls	LED Display
Heaviest Piece (when	20 kg (44 lbs) Front section

disassembled)	
Total Weight (without batteries)	44 kg (100 lbs)
Suspension	Seat Post
Battery Requirements	Two 12-volt, deep cycle (standard) Size: 2x 18 AH
Seating	Foldable, compact, vinyl-covered molded plastic
Tyres	Front: 7.62 cm x 22.9 cm (3" x 9") solid Rear: 7.62cm x 22.9 cm (3" x 9") solid

[Read more about VAT Exemption](#)

Product Gallery



