

## Sport Rider



**Brand:** Drive Medical  
**Product Code:** SPORTRID  
**Availability:** 7  
**Weight:** 117.00kg

**Price: £3,650.00**

### Short Description

This stylish Sport Rider Mobility Scooter was inspired by some of the leading designs in the Motorcycle industry.

### Description

Drive medical Sport Rider moves away from the traditional type of mobility scooter, creating desirability and excitement to the driving experience.

Using technology borrowed from the automotive industry the all round suspension provides excellent ride quality and handling, making a cruise on the Drive Medical Sport Rider Mobility Scooter both exhilarating and comfortable.

### Specifications

- Length: 65" (165cm)
- Width: 31" (78cm)
- Height: 44" (112cm)
- Weight Limit: 28 stone (180kg)
- Turning Radius: 62"
- Wheel diameter: 19.5"
- Top Speed: 8mph
- Range: 31 miles

- Total Weight: 117kg
- Battery: 100ah 12V x 2
- Suspension: Telescopic (front) and 2 x shock absorber (rear)
- Charger: 8A 24V
- Seat: Captains recline/headrest
- Motor: 1300W 4-pole
- Controller: Dynamics 160A Rhino
- Brake: EM with drum brake

**VAT RELIEF APPLIED**  
[Read more about VAT Exemption](#)

## Product Gallery



