Phone: 0117 9155253 - Email: sales@marksmobility.co.uk

Colt Executive



Brand: Pride

Product Code: COLTEX

Availability: 6 Weight: 117.80kg

Price: £4,409.00

Short Description

The Pride Colt Executive is the flagship of the Colt range of mobility scooters. A high performance 8mph road legal mobility scooter, the Colt effortlessly blends power and luxury.

Description

Colt Executive is the most stylish scooter, and comes with a wealth of high performance features as standard.

The Colt features all the high end features you could hope for: full suspension, sealed hydraulic brakes, extra large low-profile tires, advanced LED lighting, including directional signals, and quick adjust tiller.

The suspension features a mono-shock unit at the rear. This lets the Executive feel both comfortable and sporty, without the mushy feel of some large suspension scooters.

Plush captains style seat with with a pillow top that is fully adjustable and rotates for ease of access.

An excellent choice for anyone looking for comfort and style backed up by Pride's renowned reliability.

Wrap around delta tiller and LED lighting

• 25 - 30 miles range and user weight up to 28.5 stone

• Speed: 8mph

Max Weight Capacity: 28.5 stoneGround Clearance: 13.33cm/5.25"

Turning Radius: 210.cm/83"Overall Length: 142.8cm/56.25"

Overall Width: 69.2cm/27.25"

• Front and Rear Tyres: 10 cm x 33 cm (4" x 13") pneumatic

Suspension: Front and rear

• Total Weight without batteries: 117.8kg (259.lbs)

• Heaviest Pieces when disassembled (main frame): 93.21kg (205.5lbs)

• Battery Requirements: 2 x 12 Volt, deep cycle

• Battery Size: 100 AH (standard)

• Battery Weight: 27.44kg (60.5lbs) each

• Battery Charger: 8amp off board

Insurance



Product Gallery





