

Height Adjustable Kitchen Strolley Trolley with Brakes



Brand: Aidapt
Product Code: VG798WB
Availability: 15
Weight: 0.00kg

Price: £122.10

Short Description

Height Adjustable Kitchen Strolley Trolley Is a walking aid for use within the home, providing a means of transporting items from one room to another.

Description

The ideal trolley to transfer meals and other items around the home. The Strolley Trolley is height adjustable and the unique patented trays clip on and off easily for cleaning. This new Strolley Trolley has the added benefits of handle brakes for extra piece of mind for the user. The materials used are highly durable, making it a great value product for both the prescriber and the user. The castors are manufactured from polypropylene and are amazingly resistant to damage, they also incorporate a patented wheel guard to reduce hair and fibre entrapment. Supplied flat packed.

- Handy wheeled trolley is ideal for transferring meals/drinks between rooms at home.
- Suitable for use around the home
- Handle brakes for added piece of mind
- Height adjustable legs to suit the majority of users.
- Clip On/Off Trays are easy to clean and maintain

- Depth (mm): 530
- Height (mm): 840
- Width (mm): 530
- Product Dimensions (mm): 840x530x530
- Net weight (kg): 8.5
- Colour: Brown and Cream
- Configuration: With Brakes
- Footprint (WxD) Min/Max (mm): 530x530
- Height of Handle from Floor (mm): 840-1020
- Maximum Weight On Tray: 15kg
- Tray (WxD) (mm): 463x345

IMPORTANT NOTICE

This product is not weight bearing. If help is needed when walking please see our range of walking frames or for a dual purpose walking aid and trolley please see the Let's Go Indoor Rollator (M14686).

