

## Bendable Cutlery Set (3 pieces)



**Brand:** Aidapt  
**Product Code:** VM914  
**Availability:** 20  
**Weight:** 0.00kg

**Price:** £46.20

### Short Description

Dining aid set including a bendable spoon and bendable fork together with a rocker knife and utensil strap.

### Description

Bendable Cutlery Set, is a Dining Aid Set consisting of a bendable spoon and bendable fork and a rocker knife, together with a utensil strap. The bendable spoon and fork and rocker knife all have a fixed loop at the end of the wide, soft cushion handle, through which the utensil strap can be passed to secure the dining aid utensil to the user's hand.

- Set consists of a bendable spoon, bendable fork, rocker knife and a utensil strap
- Dishwasher safe up to 70 degrees
- Ideal for people with a limited grip
  
- Depth (mm): 30
- Height (mm): 210
- Width (mm): 80
- Product Dimensions (mm): 210x80x30
- Net weight (kg): 0.26
- Colour: Black Handle, Metal Head

## Product Gallery

