

Leg lifter



Brand: Aidapt
Product Code: VM949
Availability: 25
Weight: 0.90kg

Price: £11.55

Short Description

This handy personal aid enables a person to lift their legs onto a bed or into a car independently.

Description

This handy personal aid enables the leg to be lifted on or off a wheelchair or bed by pulling the handle loop. It has a wired frame inside to keep the leg lifter straight.

- Handy aid loops around your foot Allows the user to lift an immobile leg in and out of bed or into a wheelchair
- Loops around your foot and can be lifted using the handle loop
- Depth (mm): 10
- Height (mm): 830
- Width (mm): 25
- Product Dimensions (mm): 830x25x10
- Net weight (kg): 0.09
- Colour: Blue

