


Birling Bed Back Rest



Brand: Aidapt
Product Code: VG820A
Availability: 15
Weight: 2.05kg

Price: £41.24

Short Description

Birling Bed Back Rest sits users up in bed - eliminating the need for carers to lift.

Description

The Birling Bed Back Rest is an ideal and versatile adjustable back rest for resting patients or those seeking greater comfort in bed. With multiple fixed back positions giving support to the user, the Birling Bed Back Rest is covered with an easy clean fabric and comes complete with a head cushion. Supplied in a retail box.

- Can be set at 5 different angles up to 85 degrees
 - 108 kg (17 stone) weight limit
 - Easy-clean fabric and padded head cushion
 - Light but sturdy framework
 - Epoxy coated frame in white
-
- Depth (mm): 500
 - Height (mm): 280-550
 - Width (mm): 490
 - Product Dimensions (mm): 280x490x500
 - Net weight (kg): 2.05
 - Colour: Blue and white

- Depth of base: 500
- Maximum User Weight (kg): 108
- Maximum User Weight (stones): 17
- Width of Base (mm): 490

[Read more about VAT Exemption.](#)

Product Gallery

