Foam Leg Raiser

Brand: Aidapt

Product Code: VG885

Availability: 25 Weight: 0.00kg



Price: £44.55

Short Description

Foam Leg Raiser is an ergonomically-designed Foam Leg Raiser.

Description

The cotton twill covered Foam Leg Raiser has been ergonomically-designed to support swollen legs and stimulate circulation; it is ideal for use in a bed or flat surface. Constructed from a high density polyurethane foam with a generous 2 cm deep memory foam topper, the Foam Leg Raiser features a removable, washable cover.

- Removable cover for washing
- Ergonomically-designed shape
- Ideal for use in bed or flat surfaces
- · High density foam for support and comfort
- Supports swollen legs and can help to improve circulation

• Depth (mm): 170

• Height (mm): 620

• Width (mm): 415

• Product Dimensions (mm): 620x415x170

• Net weight (kg): 0.9

· Colour: White

• Depth of base: 630

• Overall Depth (mm): 630

• Overall Width (mm): 420

• Width of Base (mm): 420



Product Gallery

