

Foam Leg Raiser

Brand: Aidapt
Product Code: VG885
Availability: 25
Weight: 0.00kg



Price: £44.55

Short Description

Foam Leg Raiser is an ergonomically-designed Foam Leg Raiser.

Description

The cotton twill covered Foam Leg Raiser has been ergonomically-designed to support swollen legs and stimulate circulation; it is ideal for use in a bed or flat surface. Constructed from a high density polyurethane foam with a generous 2 cm deep memory foam topper, the Foam Leg Raiser features a removable, washable cover.

- Removable cover for washing
- Ergonomically-designed shape
- Ideal for use in bed or flat surfaces
- High density foam for support and comfort
- Supports swollen legs and can help to improve circulation

- Depth (mm): 170
- Height (mm): 620
- Width (mm): 415
- Product Dimensions (mm): 620x415x170
- Net weight (kg): 0.9
- Colour: White

- Depth of base: 630
- Overall Depth (mm): 630
- Overall Width (mm): 420
- Width of Base (mm): 420

[Read more about VAT Exemption...](#)

Product Gallery

