

## Bath Belt



**Brand:** Relaxa  
**Product Code:** CR  
**Availability:** 20  
**Weight:** 15.90kg  
**Location:** Bristol and surrounding areas

**Price: £899.00**

### Short Description

Bath lifts fit inside the bath and can lift or lower the bather from near the bottom of the bath up to the height of the bath rim.

### Description

The bath belt provides independence for those that find the bath difficult to enter or leave.

Comprises of an inconspicuous wall-unit and a discreet bracket at floor level.

With the easy to use hand set, a strong broad belt feeds from the wall-unit. This is easily secured to the floor bracket. The belt is then made taut with the touch of another button to create a secure, very stable seat.

The easy to use hand set then slackens the belt, lowering you gently down to enjoy a totally normal bath, and to come back up is every bit as straight forward.

### Specifications

- The unit has a thermometer to test the temperature of the water
- Two 12 volt rechargeable lead acid batteries provide a 24 volt DC supply to the motor
- Indicator lights in green and red shows the battery is in operation and the

battery condition

- A managed 24 volt 300MA charger is provided
- Maximum user weight 127 kgm (20 stone)
- "UP" and "DOWN" buttons are located on the front of the bathlift. A sealed hand control with "UP" and "DOWN" buttons is also provided
- An intergral handle is provided but this can be removed or re-sited as required
- Made from Re-inforced Polypropylene

## Dimensions

The unit is slim and fits neatly at the side of the bath.

- Length 55.5cm
- Height 32.5cm
- Depth 11.8cm

## Weight

The unit itself weighs 11kgms. The box when packed to include other items.

- Weight 15.9kg
- Height 32.5cm
- Depth 11.8cm

- The complete unit is manufactured in the UK.
- Maximum user weight up to 20 stones in weight.

[Read more about VAT Exemption...](#)

## Product Gallery

