

Folding Ramps



Brand: Drive Medical
Product Code: R009 6FT / ROO9 3FT
Availability: 15
Weight: 0.75kg

Price: £119.00

Short Description

Folding Ramps to aid access up and down steps and kerbs.

Description

The Folding Ramps are designed to aid access up and down steps and kerbs.

- High quality user friendly alloy ramp
- Provides a perfect and safe solution for access
- Ideal for Wheelchairs, Powerchairs and Scooters
- Features raised edges for added safety

Product Features

- They provide a perfect and safe solution for disability and mobility access.
- Ideal for Wheelchairs, Powerchairs and Scooters
- The two hinged tracks simply unfold and can be positioned accordingly for maximum safety and support with an angled lip at the top of the ramp for added security.
- Ideal for use by retail shops on the high street in accordance with the DDA (Disability Discrimination Act).
- For temporary or permanent installation - two holes included for

permanent installation.

- These folding aluminium ramps are versatile, lightweight and easy to use.
- Features raised edges for added safety and peace of mind.
- Constructed from lightweight aluminium and includes a central plastic handle for easy carrying and transportation.
- Two aluminium sections 35cm and 37cm wide to fold to a compact size.
- High quality user friendly alloy ramp.
- A wide anti-slip grip surface ensures safe use in most weather conditions.
- Durable design to support both wheelchair/mobility scooters and the occupant with a weight capacity of up to 272kg SWL (safe working load).

R009-3FT Folding Ramp:

- LENGTH 3ft (0.9m)
- WIDTH 72cm (28")
- LENGTH FOLDED 91cm (36")
- WIDTH WHEN FOLDED 32cm (12.5")
- WEIGHT 9kg (20lb)
- WEIGHT CAPACITY 272kg (43st)

R009-6FT Folding Ramp:

- LENGTH 6ft (1.8m)
- WIDTH 72cm (28")
- LENGTH FOLDED 183cm (72")
- WIDTH WHEN FOLDED 32cm (12.5")
- WEIGHT 18kg (40lb)
- WEIGHT CAPACITY 272kg (43st)



Product Gallery

