

Bath Lift Bellavita



Brand: Drive Medical
Product Code: 470100252
Availability: 10
Weight: 9.30kg

Price: £499.99

Short Description

The reclining Bath Lift Bellavita, with ultra modern design featuring a unique front aperture for improved personal hygiene.

Description

This ultra modern bathlift takes the Bellavita away from an austere and clinical look and is more in keeping with today's requirements. Developed with the help of Nursing Professionals, it includes several excellent new features.

- The lightest bath lift on the market at only 9.3kg
- At 6cm the Bellavita is the lowest bathlift available
- Backrest reclines to a market leading 50°
- Reaches a record seat height of 48cm
- Lifetime guarantee on frame & motor
- Ergonomically designed for greater comfort and safety for the user
- Easy to install and operate
- User friendly charger
- Quick charge lithium battery technology
- Built in safety controls
- Lightweight, fully waterproof floating hand control
- Washable hygienic covers available in white, blue or grey
- Hygienic cutaway and slimline profile for easy handling

- Easy to dismantle for storage and transportation
- 4 quick release suckers with an option for an additional two for an anti-slip bath base
- Standard Bellavita weight capacity - 22st (140kg)
- Bellavita HD weight capacity - 25st (160kg)
- 2 yrs guarantee on battery/handset
- Headrest available as an optional extra

Code	Length	Width	Seat Width	Seat Depth	Seat Length	Backrest Height	Baseplate Length	Baseplate Width
470100	83cm (33")	70cm (28")	35cm (14")	42cm (17")	57cm (22")	64cm (25")	56cm (22")	31cm (12")

[Read more about VAT Exemption](#)