

Light weight 3 wheeled



Brand: Drive Medical
Product Code: TW008
Availability: 30
Weight: 0.00kg

Price: £119.00

Short Description

The Walking Frame Apron Bag is a lightweight and hugely useful item that can easily be folded and slipped into a bag or pocket when not in use. Walking aid developed to keep its weight to an absolute minimum providing the features of an outdoor Triwalker.

Description

Light weight 3 wheeled walking aid was developed to keep its weight to an absolute minimum whilst still providing the features and benefits associated with an outdoor Triwalker of this class.

Specifications

- Lightweight solid 7" wheels for indoor or outdoor use
- Comes standard with vinyl bag
- Arthritic friendly loop lock brakes
- Handles are adjustable in height
- Easy, one hand folding
- Dimensions (cm): 61W x 56D x 83 – 93H (24" x 25.5" x 33"–37")
- Total Weight 4.7kg (10lb)
- Weight capacity 115kg (18 stone)

VAT RELIEF APPLIED
[Read more about VAT Exemption.](#)