

## Steel triwalker with bag, basket, tray



**Brand:** Drive Medical  
**Product Code:** TW015  
**Availability:** 30  
**Weight:** 0.00kg

**Price: £115.50**

### Short Description

The triwalker with bag, basket and tray features a durable yet lightweight steel frame making it easy to manoeuvre and transport.

### Description

This triwalker with bag, basket and tray features a durable yet lightweight steel frame making it easy to manoeuvre and transport. The arthritic friendly looped cable brakes can be easily activated by the user.

- Sturdy steel frame construction is durable
- Looped cable brakes are adjustable
- Handle height are adjustable
- Easy fold system for storage and transportation

### Product Features:

- Sturdy steel frame construction is durable yet lightweight for manoeuvrability and transportation
- Arthritic friendly, looped cable brakes are adjustable, easily activated and lockable for safety
- Handle height adjusts between 84cm – 99cm (33" - 39" ) to suit user requirements

- 19cm (7.5") lightweight wheels and arthritic friendly cable brakes
- Easy fold system for storage and transportation
- Includes a large basket and vinyl storage bag for carrying shopping and other essentials
- Includes a practical tray with raised edges to prevent dropping items

<b>Specifications:</b> WIDTH	62cm (25")
DEPTH	60cm (24")
HEIGHT MIN	84cm (33")
HEIGHT MAX	99cm (39")
WEIGHT	7kg (15lb)
WEIGHT LIMIT	124KG (20ST)

[Read more about VAT Exemption](#)

## Product Gallery

