

Light weight walking aid with seat and basket



Brand: Drive Medical
Product Code: R6/8
Availability: 10
Weight: 0.00kg

Price: £129.99

Short Description

Walking aid developed to keep its weight to an absolute minimum providing the features of an outdoor Tri-Walker

Description

Light weight 4 wheeled walking aid with seat and basket is light to carry and light on your pocket, this New Aluminium Rollator offers many upmarket features at a price you can afford.

Specifications

- Very lightweight at less than 7kg
- Height adjustable, 33"- 38"
- Anatomic handgrips
- Arthritic-friendly loop cable brakes
- Padded flip-up seat
- Under seat storage bag
- Easy-to-use folding mechanism
- 7"(18cm) diameter solid tyres (R8)
- 6" (15cm) diameter solid tyres (R6)
- Curved flip-up backrest for comfort and easy storage
- Choice of two metallic colours; Vivid Blue or Spirit Red.

Product Gallery

