

Walking Stick Ice Grip



Brand: Drive Medical
Product Code: WSA006RT
Availability: 30
Weight: 0.00kg

Price: £9.99

Short Description

Walking Stick Ice Grip provides grip in snow or icy weather conditions.

Description

Specifications:

- Five pronged Ice Grip assists users when using their walking aids in snow or icy weather conditions
- Suitable for most walking sticks and crutches with tubing / ferrule sizes between 19mm and 22mm
- Ice grip can be 'flipped up' when not in use into a locked position to allow walking stick to be used as normal
- Supplied with 2 pairs of screws; shorter screws to install the ice grip on narrow tubing and longer screws for wider tubing

 [Read more about VAT Exemption](#)

Product Gallery

