Phone: 0117 9155253 - Email: sales@marksmobility.co.uk

Pvc Handledelbow Crutches



Brand: Aidapt

Product Code: VP147/8

Availability: 30 Weight: 0.00kg

Price: £39.99

Short Description

This lightweight aluminium elbow crutch has a wedge-shaped PVC handle for improved grip, along with a non-slip rubber foot.

Description

Pvc Handledelbow Crutches also features double height adjustment for comfort and flexibility. Priced per pair of elbow crutches, please ensure that you select the correct size from the available options when placing your order..

- Lightweight aluminium elbow crutch ideal for use during rehabilitation
- Double height adjustment with 14 plus 3 height settings
- Soft touch handgrip for greater comfort
- Maximum user weight limit 160kg (25 stone)
- Rubber foot to prevent slips

Specifications

Code	VP147	VP148
Colour	Silver/Grey	Silver/Grey
Depth (mm	155	155
Height (mm)	775-1055	870-1155
Height to Handgrip (mm)	575-815	680-925

Inverter	No	No
Length (Collapsed) (mm)	775	870
Maximum User Weight	160	160
(kg)		
Maximum User Weight	25	25
(stones)		
Net weight (kg)	1.4	1.4
Number of Height Settings	s 17	17
Overall Height (mm)	1055	1155
Product Dimensions (mm)	n/a	680x100x155
Size	Medium	Adult (Large)
Width (mm)	100	100
Read more about VAT Ex	emption	