

Travelite 18 stone weight limit wheelchair



Brand: Drive Medical
Product Code: Tc005
Availability: 8
Weight: 8.00kg

Price: £199.00

Short Description

The new TraveLite transit wheelchair from Drive Medical is totally portable, lightweight and comes with a strong carry holdall.

Description

Travelite 18 stone weight limit wheelchair is totally portable, lightweight and comes with a strong carry holdall.

Features

- Aluminium frame in metallic gun metal grey
- Half folding back mechanism
- 18" square nylon seat canvas with lift up straps
- Parking brakes
- Quadruple steel cross brace
- Padded, flip up desk arms
- Height adjustable footrests
- 8" solid tyres

Specifications

- Width open: 57cm

- Width closed: 29cm
- Length: 84cm
- Height: 95cm
- Seat Width: 45cm
- Seat Length: 45cm
- Seat Height: 48cm
- Weight: 9.3kg
- User Weight Capacity: 18 stone

VAT RELIEF
[Read more about VAT Exemption...](#)

Product Gallery

