# Phone: 0117 9155253 - Email: sales@marksmobility.co.uk

## Travelite 18 stone weight limit wheelchair



**Brand:** Drive Medical **Product Code:** Tc005

Availability: 8 Weight: 8.00kg

Price: £199.00

#### **Short Description**

The new TraveLite transit wheelchair from Drive Medical is totally portable, lightweight and comes with a strong carry holdall.

#### **Description**

Travelite 18 stone weight limit wheelchair is totally portable, lightweight and comes with a strong carry holdall.

#### **Features**

- Aluminium frame in metallic gun metal grey
- Half folding back machanism
- 18" square nylon seat canvas with lift up straps
- Parking brakes
- Quadruple steel cross brace
- Padded, flip up desk arms
- Height adjustable footrests
- 8" solid tyres

### **Specifications**

Width open: 57cm

• Width closed: 29cm

Length: 84cmHeight: 95cm

Seat Width: 45cmSeat Length: 45cmSeat Height: 48cm

• Weight: 9.3kg

• User Weight Capacity: 18 stone



**Product Gallery** 



