

Aktiv X2 Lite Aluminium



Brand: Aktiv Wheelchairs
Product Code: X2
Availability: 8
Weight: 8.00kg

Price: £289.00

Short Description

The Aktiv X2 is a lightweight durable wheelchair for occasional use.

Description

The Aktiv X2 is a lightweight durable wheelchair. It has a lightweight aluminium folding frame and is easy to lift, store and transport. The armrests flip back and can be removed easily as standard and have comfortable polyurethane armrest pads. The backrest folds down and the leg rests and rear wheels are removable by easy quick release mechanisms; therefore the chair can be quickly and easily disassembled for easy storage in a car boot.

Features

- Folding Aluminum frame
- Lightweight aluminium wheelchair
- Available in Transit (12" rear wheel) or Self Propel (24" rear wheel) models
- Dual axle allows for easy conversion between transit & self propel models
- Flip back, detachable armrest with PU arm pad
- Plastic side panel clothes guard
- Swing away leg rests with black plastic foot plate & heel loops
- Black nylon design seat and back upholstery

- 7" Grey PU front castors
- Solid Tyres
- Press button release safety belt

Specifications

- Seat Width: 41cm (16"), 46cm (18"), 51cm (20")
- Seat Depth: 41cm (16")
- Seat Height: 48cm (19") from ground the seat canvas)
- Backrest Height: 45cm (18")

Weight Specs

- Total Weight: 13.2kg
- Transport Weight: 9.1kgs
- Max Weight: 113kg / 17.75st

This chair has passed the crash test ISO7176-19:2001 for occupied transportation

*A Crash Tested Wheelchair has a frame which is manufactured to comply with test ISO7176-19:2001 for occupied transportation



Product Gallery

