

## Norfolk chemical commode



**Brand:** Aidapt  
**Product Code:** VR235  
**Availability:** 15  
**Weight:** 0.00kg

**Price: £79.20**

### Short Description

The stackable Norfolk Commode Chair can accommodate 5, 10 and 20 litre commode pails...

### Description

The stackable Norfolk Commode Chair is the latest addition to our range of commodes and is ideal as a bathroom commode or bedroom commode chair. It is height adjustable and has the added benefit of an adjustable sliding pan rack, this enables the Norfolk Commode Chair to accommodate 5, 10 and 20 litre commode pails. As with all our commodes, the seat and back are vinyl welded for hygiene and safety.

- Ideal for users who struggle to get to the bathroom
- Height adjustable to suit the majority of users
- Fixed padded back and detachable padded overseat for comfort
- 5, 10 and 20 litre pails available
- Supplied with a 5 litre pail with lid as standard
- 154 kg (25 stone) weight limit
- Capacity (litres): 5
- Colour: Black and brown
- Height (mm): 800-900
- Height from Floor to Top of Back (mm): 800-900

- Height from Floor to Top of Seat (mm): 480-570
- Height of Arms from Seat (mm): 225
- Maximum User Weight (kg): 154
- Maximum User Weight (stones): 25
- Seat Depth (mm): 350
- Width between Arms (mm): 505

 [Read more about VAT Exemption...](#)